


FIG 1

Lightweight 3 Wheel rollator with adjustable handles and a basket.

Technical Specification

- Maximum user weight limit: 136kg/300lbs
- Maximum basket weight load: 5Kg
- Maximum bag weight load: 5Kg

ASSEMBLY INSTRUCTIONS:



WARNING

After locking the knob ensure that left and right handles are the same height.




Brake



Lock the brakes



Release the brakes



Brake lock status



USER INSTRUCTIONS:

- Ensure the Tri Wheel Rollator is assembled in as shown in Fig 1 above.
- Stand upright behind the rollator with your shoulders relaxed and your arms by your side, looking straight ahead.
- In this position, the handgrips of the rollator should be adjusted to a height level with your wrists. This will position your arms at an approximate 20° to 30° bend when using the rollator.
- Check that the holes in the handle bars and the frame are aligned.
- Only use the rollator as a walking support.
- Do not use the rollator on stairs or an escalator.
- Exercise caution when there is a heavy load in the bag (max. 5Kg) or basket /tray (max.5Kg)
- Loads should only be transported using the bag or basket/tray (where fitted). Do not hang heavy items from the handles as this may affect the stability of the Rollator.
- Do not use the rollator to walk backwards.
- Caution must be taken when negotiating curbs and other obstacles, or when using on sloping, uneven or slippery surfaces. Failure to do so may result in a serious risk of a fall or injury.
- Do not perform any adjustments to the rollator while it is in use.
- Wheels must be in contact with ground / floor at all times.
- When using the rollator in a stationary position, the brakes must be locked.

CLEANING INSTRUCTIONS:

- Clean and disinfect regularly
- Store in a dry and ventilated environment.

MAINTENANCE AND SAFETY:



- Always check whether the brake is flexible and whether the hand grip, folded part, is fully open in place.
- Check whether the wheels are turning and are not worn or shaking.
- Ensure the wheels and bearings are properly maintained.
- Check whether anti-slip rubber sleeve is damaged or not, if it is worn, please replace it.
- Regularly check for wear and tear and that all screws and knobs are securely fastened.
- Check brakes are working satisfactorily by squeezing the brakes and by braking to a stop.

DISPOSAL:

- Dispose of all packaging, paper, cartons and plastic, plastic bags in accordance with your local recycling regulations.
- At the end of the product's lifespan please dispose of it at an authorised household waste recycling centre.

