

Air Compression Leg Massager

K632 v01: 18/03/22

Please read and retain these instructions for future reference

vivadia



Dimensions

(H) 40 x (L) 28cm

Length of mains lead: 170cm

- AIR COMPRESSION MASSAGER THAT PREVENTS AS WELL AS EASES – HEAT IMPROVES CIRCULATION, MASSAGE PROVIDES RELIEF
- USE ONCE OR TWICE A DAY FOR 10 MINUTES EACH TIME TO RELIEVE LEG FATIGUE. MULTIPLE MASSAGE MODES BRING YOU A DIFFERENT MASSAGE EXPERIENCE.
- AIR COMPRESSION CHAMBERS CONTRACT AND MASSAGE MUSCLES EFFECTIVELY RELIEVING SWELLING, LEG FATIGUE, AND BRINGING INCREDIBLE COMFORT TO THE LEGS.
- FOR ANKLES, CALVES AND FEET A LOWER LEG MASSAGE
- HOOK & LOOP DESIGN SO YOU CAN GET A COMFORTABLE FIT
- MAX CIRCUMFERENCE: 50cm

Technical details

Input Voltage: 100-240Vac, 50Hz-60Hz, 1.0A Max

Output Voltage: 12.0V=2.0 A, 24W

Default time:20min.

Please consult your GP if you are pregnant, suffer from skin infections, high blood pressure, use a pacemaker, or another serious medical condition.

Safety Precautions

- DO NOT USE VOLTAGE BEYOND THE REQUIREMENT OF THIS PRODUCT.
- DO NOT ATTEMPT TO REPAIR OR DISMANTLE THE PRODUCT.THERE ARE NO USER-SERVICEABLE PARTS.
- DO NOT USE THE PRODUCT IN A PLACE WHERE IT MAY BECOME WET OR DAMP. IF THE UNIT BECOMES WET, DO NOT USE.
- STOP USING THE PRODUCT IMMEDIATELY IF YOU FEEL UNCOMFORTABLE OR ANY ABNORMALITIES OCCUR.
- DO NOT USE THIS PRODUCT IN A HIGH-TEMPERATURE ENVIRONMENT
- DO NOT USE THIS PRODUCT IF YOU ARE CURRENTLY PREGNANT OR MENSTRUATING.
- DO NOT PLACE THE UNIT IN A PLACE WHERE IT MAY BE IN CONTACT WITH FLAMMABLE ITEMS, EXPLOSIVE ITEMS. OR CHEMICALS.
- DO NOT TOUCH THE UNIT WITH ANY SHARP OBJECT OR RISK DAMAGE.
- DO NOT USE THIS PRODUCT IF ANY OF THE BELOW APPLY TO YOU:
 - SKIN IS INFECTED, SWOLLEN, BURNT,OR OTHERWISE DAMAGED



**The Enterprise
Department**

UK Head Office, 11 Bridge Street, Bishop's Stortford, CM23 2JU

Tel: 0844 4824400*, Fax: 01279 756595

www.coopersofstortford.co.uk

EU address: D.S.B. LD, 160 Bd. de Fourmies, 59100 Roubaix, France


Tel: (+44) 1279 701269

*UK Calls cost around 7 pence per minute plus your telephone company's network access.
Please check with your telephone supplier.


- YOU HAVE A PACEMAKER OR OTHER MEDICAL IMPLANT.
- YOU ARE SUFFERING FROM HEART DISEASE.
- YOUR BONES ARE FRAGILE OR INJURED.
- YOU ARE SUFFERING FROM CANCER.
- YOU ARE RECEIVING TREATMENT OR YOUR BLOOD PRESSURE IS NOT NORMAL.
- YOU ARE SUFFERING FROM OSTEOPOROSIS.

Controls



1. Press the  key to enter the slow horizontal mode automatically, the air pressure is in the low position by default, and the heating function is in the off state.
2. Press the mode switch to cycle through the three modes (relieve, relax, and smart).
3. Press the air pressure switch, and the three strengths (high, middle, low) are switched in cycles.
4. Press the heating switch to switch the temperature in 2 levels (high level and low level) cyclically.



5. Press the  again to shut down and end the work.



Operation

1. Slip on each massager, keeping your feet flat on the floor.
Note that one is marked L and the other R for left and right foot.
2. Use the strap to tightly secure your foot into the massager.
3. Use the strap to tightly secure your leg into the massager.
4. Move your foot to make sure that the massager is securely fastened.
5. Insert the two air pressure hoses into the controller.
6. Plug the adaptor into a wall socket.
7. Sit comfortably. Do not cross your legs.
8. Press the power button on the control panel. When switched on the LED will light up.
9. Adjust the air pressure strength and mode which will be shown on the LCD screen.



NOTE: the device will automatically shut off after 20 minutes
For best results, do not use the leg massager for more than an hour.

Care and maintenance

- Store the product in a dry place, keep away from moisture and direct sunlight.
- Wipe clean with a dry towel.

Troubleshooting

Problem	Cause	Solution
The device does not work.	The unit is not connected to a power source.	Make sure the unit is fully connected to the power.
	The power button was not pressed.	Press the power button
The device suddenly stopped working.	The unit was running more than 15 minutes.	Restart the unit to resume use

PLUG WIRING (UK & IRELAND)

This appliance is supplied with a mains adaptor/charger with BS1363 plug pins.


DO NOT use If the socket outlets in your home or office are not suitable for this product's mains adaptor/charger plug pins. Attempts to insert the mains adaptor/charger into the wrong socket is likely to cause electric hazard or damage.

DO NOT use with any other mains transformer/charger than the one that is supplied with this product.

DISPOSAL

- Coopers of Stortford use recyclable or recycled packaging where possible.
- Please dispose of all packaging, paper, cartons, packaging in accordance with your local recycling regulations.
- Plastics, polybags – Contains the following recyclable plastic.



Code & Symbol	 OTHER
Type of Plastic	ABS, TPU, Silicone, Polyester, Nylon
Commonly used for	CD's, crisp packets, various flexible packaging, baby bottles, sunglass lenses
Notes	Not often or widely recycled

- At the end of the product's lifespan please check with your local council authorised household waste recycling centre for disposal.

UK
CA

CE



The Enterprise Department

UK Head Office, 11 Bridge Street, Bishop's Stortford, CM23 2JU

Tel: 0844 4824400*, Fax: 01279 756595

www.coopersofstortford.co.uk

EU address: D.S.B. LD, 160 Bd. de Fourmies, 59100 Roubaix, France

Tel: (+44) 1279 701269

*UK Calls cost around 7 pence per minute plus your telephone company's network access.
Please check with your telephone supplier.