

# Folding Dog Steps

G407

# Instructions for Use

Please read and retain these instructions for future reference

v001: 08/11/19

## OPEN



**Max weight load:**  
50kg


**Step Tread:**  
35cm x18cm

## FOLDED



**Dimensions when folded:**  
46cm x41cm x12cm  
Open height adjustable from 50cm to 70cm

PERFECT FOR OLDER DOGS WHO MAY HAVE ARTHRITIS OR JOINT PROBLEMS  
PORTABLE AND COLLAPSES FLAT SO THEY'RE EASY TO TRANSPORT  
HANDY CARRY BAG FOR EASY STORAGE  
ANTI-SKID STEPS  
LIGHTWEIGHT YET STRONG FRAME  
WILL FIT SMALL AND LARGER SIZED CARS  
ALSO CAN BE USED INDOORS IF YOUR DOG NEEDS A HELPING HAND ONTO A BED OR SOFA



\*Calls cost 7 pence per minute plus your telephone company's network access charge

Head Office, 11 Bridge Street,  
Bishop's Stortford CM23 2JU  
Tel: 0844 482 4400\*  
Fax: 01279 756 595  
[www.coopersofstortford.co.uk](http://www.coopersofstortford.co.uk)

**Coopers**  
OF STORTFORD



### **SAFETY WARNING**



The folding pet step is designed to support a maximum static capacity of 50kg. Exceeding this capacity may cause damage to the product or cause injury to the pet.  
Beware of pinching hazards when setting up or packing up the folding dog steps. Keep fingers away from pivoting frame.

The folding dog steps are not designed to be used by people.  
Keep Folding Dog Steps away from any naked flame or heating device.  
Always remember to leave animals in well ventilated areas with available drinking water.

#### **DISPOSAL**

- Dispose of all packaging, paper, cartons, plastic and plastic bags in accordance with your local recycling regulations.
- At the end of the product's lifespan please dispose of it at an authorised household waste recycling centre.